
**BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**

## ◆ DEVICE DESCRIPTION :

The LAGIS Cap for wound retractor is a single use device sterilized by ethylene oxide and must be used with the compatible size of LAGIS Wound Retractor. It comes with a needle to facilitate puncturing on the surface and to allow easy insertion of trocars.

## STRUCTURE

## ◆ INDICATIONS/INTENDED PURPOSE :

- The LAGIS Cap for wound retractor must be used with the compatible LAGIS Wound Retractor, whose combination could be taken as a single-port platform for inserting the trocars while maintaining pneumoperitoneum during the endoscopic surgery.
- Expected clinical benefits: provide a flexibly single-port platform for the number and position of trocar placement during laparoscopic surgery and can be rapidly converse from laparoscopic surgery to open surgery, if necessary.
- Patient target group: The LAGIS Cap for wound retractor is intended for adults who need endoscopic surgery.

## ◆ CONTRAINDICATIONS :

The LAGIS Cap for wound retractor is not intended for use when minimally invasive techniques are contraindicated.

## ◆ INSTRUCTIONS FOR USE :

The LAGIS Cap for wound retractor shall be used in combination with the following compatible models of the LAGIS Wound Retractor:

| LAGIS Cap for wound retractor for Wound Retractor |                                  |
|---------------------------------------------------|----------------------------------|
| Model                                             | Compatible LAGIS Wound Retractor |
| CP-60ES                                           | WR-60ES, WR-60ES-H               |
| CP-80S                                            | WR-80S                           |
| CP-110M                                           | WR-110M-A                        |

Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised.

### Installation and Assembly:

- Installation of trocars: Puncture a hole at the puncturing area with the needle to facilitate insertion of a trocar by rotation motion and continuous controlling downward pressure. The number of holes depends on the required quantity of trocars during the surgery.
- Assembly: Once the LAGIS Wound Retractor is securely installed, attach the LAGIS Cap for wound retractor with inserted trocars to the attachment ring and make sure it is fully embedded in the groove of the cap. After assembly, adjust the depth of inserted trocars to prevent the trocars from damaging the organs or the retraction sleeve of the wound retractor.
- Insufflation: Connect the insufflation tubing to the stopcock located at the side of the trocar cannula to establish pneumoperitoneum. If establishing pneumoperitoneum is failed, check if the LAGIS Cap for wound retractor is tightly attached to the attachment ring.

- To insert an additional 5 mm trocar during the surgery, perform step (1) and avoid damaging the organs or the retraction sleeve of the wound retractor while inserting.  
To insert an additional 10/12 mm trocar during the surgery, first of all, perform desufflation and detach the LAGIS Cap for wound retractor from the attachment ring. Afterwards, remove the inserted 5 mm trocar and replace with the 10/12 mm trocar by following step (1) and step (2).
- Disassembly: Once the surgery is completed, withdraw the inserted trocars and detach the LAGIS Cap for wound retractor from the attachment ring carefully.

## ◆ WARNINGS AND PRECAUTIONS :

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- This device is designed to be used with the LAGIS Wound Retractor only.
- If trocars or accessories from different manufacturers are used together in a procedure, verify compatibility prior to performing the procedure.
- The 10/12 mm trocars are not allowed to be inserted directly into the LAGIS Cap for wound retractor, or they may compromise the integrity of the device and lead to air leakage.
- Do not insert the trocars too closely, or it may lead to air leakage.
- Do not withdraw the inserted trocars during the procedure, or the air may leak out from the puncturing site and lose the pneumoperitoneum.
- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or reesterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in this instructions for use.
- Improper assembly or operation of the device may result in tissue or organ injuries from mechanical damage to the patient or operator.

## ◆ ENVIRONMENTAL CONDITIONS FOR STORAGE :

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

## ◆ EXPLANATION OF SYMBOLS :

|  |                     |  |                                                                                             |  |                                                                    |
|--|---------------------|--|---------------------------------------------------------------------------------------------|--|--------------------------------------------------------------------|
|  | Temperature limit   |  | Caution                                                                                     |  | Authorized representative in the European Community/European Union |
|  | Do not re-sterilize |  | Manufacturer                                                                                |  | Not made with natural rubber latex                                 |
|  | Do not re-use       |  | Consult instructions for use or consult electronic instructions for use                     |  | Do not use if package is damaged and consult instructions for use  |
|  | Date of manufacture |  | Single sterile barrier system                                                               |  | Sterilized using ethylene oxide                                    |
|  | Medical device      |  | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. |  |                                                                    |

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