**BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.**

## ◆ DEVICE DESCRIPTION:

The LAGIS Irrigation Pump is intended to pressurize an irrigation solution bag to support high pressure or high flow of irrigation solutions during surgery. It utilizes pressurized gas to provide irrigation pressure up to 700 mmHg. The user controls the flow rate using a roller clamp on the irrigation set. The maximum flow rate is determined by the viscosity of the fluid, the characteristics of the irrigation set, and the pressure applied.

## PRODUCT MODEL

IP-1000

IP-2000

## ◆ INDICATIONS/INTENDED PURPOSE:


- The LAGIS Irrigation Pump has application in general surgical procedures for rapid delivery of irrigation solutions where flow rate precision is not required.
- Expected clinical benefits: The device applied to clear debris or blood can be helpful to clearly identify the origin site of bleeding during an operative procedure. It may also reduce the rate of surgical site infections through the removal of dead or damaged tissue, metabolic waste, and wound exudate.
- Patient target group: the LAGIS Irrigation Pump is intended for adults under the surgery that requires a rapid administration of fluid products.

## ◆ CONTRAINDICATIONS:

The device is not intended for infusion of blood, IVs, cystoscopy, or for hysteroscopy.

## ◆ ILLUSTRATION OF IRRIGATION PUMP COMPONENTS

- Control Switch
- Pressure Regulator Knob
- Rapid Deflation Switch
- Latch
- Cover
- Air Source Connector
- Air Bladder (behind Cover)
- Bag Hanger
- Spring
- Pressure Gauge

## ◆ INSTRUCTIONS FOR USE:

### SPECIFICATION OF THE LAGIS IRRIGATION PUMP

|                       |                                                                                                                                   |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Input                 | Compressed air, nitrogen (Do not exceed 60 psi.)                                                                                  |
| Gauge accuracy        | ± 30 mmHg                                                                                                                         |
| Operating Pressure    | 0-700 mmHg                                                                                                                        |
| Compatible Fluid Bags | Flexible bags, supplied by manufacturers are compatible with Irrigation Pump. Rigid and semi-rigid containers are not compatible. |

## OPERATING INSTRUCTIONS

- Fix Irrigation Pump to IV pole outside the sterile field.
- Flip Control Switch to "DEPRESSURIZE" before opening the Cover.
- Connect the Air Source Connector to the compressed nitrogen or air bottles/tanks for applying pressure to the Irrigation Pump.  
Note: Connecting to a source other than those recommended, including the fluid warming system, may result in incorrect operating pressures or damage to internal components. Do not connect to nitrous oxide or other anesthetic gases. There is no indication on the Irrigation Pump that a gas pressure source has been connected.
- Open the Cover and hang the fluid bag on the Bag Hanger.
- Spike the fluid bag with an irrigation set. Use the irrigation set per the instructions for use of the irrigation set.
- Close the Cover.  
Note: If the air bladder is partly inflated, turn on the rapid deflation switch to deflate the air bladder. It will exhaust residual air from the air bladder so the cover can be closed and latched.

- Secure the Latch.

Note: Ensure the cover is latched before pressurizing.

- Flip the Control Switch to "PRESSURIZE" and adjust the pressure by pulling the Pressure Regulator Knob and rotating clockwise to increase, and counter-clockwise to decrease. Push the Pressure Regulator Knob to lock the pressure.

Note1: Do not operate the LAGIS Irrigation Pump at a pressure of 700 mmHg. The pressure has been set at the factory to approximately 300 mmHg.

Note2: When adjusting the pressure, turn the Pressure Regulator Knob clockwise to increase pressure and counterclockwise to decrease pressure. When decreasing pressure, the pressure shall be adjusted to slightly lower than the target pressure and then increased to the target pressure so that the target pressure can be maintained in actual use.

- Adjust the roller clamp on the irrigation set to achieve the desired flow rate.

Note: The LAGIS Irrigation Pump controls the pressure applied to the fluid bag. There are no alarms provided for empty bags, occlusions, or excessive flow. It is the user's responsibility to ensure proper fluid flow, including monitoring for excessive flow or occlusion.

- Change the fluid bags: Flip Control Switch to "DEPRESSURIZE" and Rapid Deflation Switch to "ON". After the Air Bladder is deflated, open the Cover and replace the fluid bags.

## CLEANING

The exterior of the pump should be wiped down with a damp cloth and germicide. Failure of the door and other parts may be caused by the application of cleaning solutions containing alcohol or solvents. No cleaning is required inside the pump. The use of other cleaning agents may result in premature failure of the product. Do not immerse the Irrigation Pump in liquids or use alcohol or solvents to clean. Do not autoclave.

## PREVENTIVE MAINTENANCE

- Monthly, the user should:
  - Examine the surface for cracks that may occur (e.g. due to the use of improper cleaning solutions). If cracks are observed, remove the Irrigation Pump from service and replace the damaged part(s).
  - Check gauge for zero reading (within one division) with no pressure applied. A non-zero reading is an indication that Irrigation Pump may have been dropped or otherwise mishandled and the unit should be thoroughly checked.
  - Check the regulator operating range of 0 to 700 mmHg. If the full pressure cannot be achieved, check for air leaks in tubing or air bladder.
- Annually, in addition to monthly checks, it is recommended to users:
  - Verify maximum pressure is between 600 and 700 mmHg with 60 psi input.
  - Check tubing and connections for leaks and signs of wear.

Note: The calibration cannot be performed by users. Factory calibration and replacement components are available if required.
- The user should replace the cover and air bladder if any cracks or damages are observed.

## DISPOSAL

The primary plastic materials used in the Irrigation Pump are Polycarbonate and Nylon. The method of disposal of the Irrigation Pump at the end of life should be appropriate for these materials according to the relevant national laws.

## ◆ WARNINGS AND PRECAUTIONS:

- The storage period shown on the label is the period of the entire equipment. The air bladder should maintain its functionality for use 10000 times, and the air bladder will not deflate when the times limit is reached. If there is a problem with the air bladder or other components, please refer to the TROUBLESHOOTING section or contact the manufacturer for the service required.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Improper assembly or operation of the device may result in prolonged surgical time and undesirable tissue injuries to the patient.
- To avoid the patient or user hurt, do not use the device when it is damaged or there is any component loosening.
- The use of excess irrigation fluid and irrigation time can lead to electrolyte imbalance, which is one of the complications, as this may cause prolonged surgical time and unintended hurt to the patient.
- High-pressure irrigation also used for hydro-dissection under an operative procedure may result in the patient hurt.
- If air sources or accessories from different manufacturers are used together with the LAGIS Irrigation Pump, verify compatibility prior to performing the procedure.
- In the event of malfunction of the device, changes in performance that may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in the instructions for use.
- The device are provided non-sterile, reusable and are to be used outside the sterile field only. The device must be thoroughly cleaned before the first use and after every use.
- When adjusting the pressure, it is a normal phenomenon that peculiar sound is produced by the Pressure Regulator Knob.
- Do not lift the Pressure Regulator Knob up when adjusting to the lowest pressure. Otherwise, this will cause idling. This phenomenon does not occur under normal use.
- When Pressure Regulator Knob is depressed, it will cause the set pressure become incorrect. The pressure can be increased or decreased in advance, which is to achieve the final one is as the expected value of pressure.

## TROUBLESHOOTING

| Problem                       | Cause                                     | Solution                                                             |
|-------------------------------|-------------------------------------------|----------------------------------------------------------------------|
| No Irrigation                 | Pump not hooked up to gas source          | Connect pump to gas source                                           |
|                               | Unit is not "ON"                          | Flip control switch to "PRESSURIZE"                                  |
|                               | Irrigation pressure is set to zero        | Pull knob and turn clockwise to reach desired pressure. Push to lock |
|                               | Rapid deflation switch is "ON"            | Flip rapid deflation switch to "OFF"                                 |
|                               | Kinked tubing                             | Straighten or replace tubing                                         |
|                               | Fluid bag is empty                        | Replace fluid bag                                                    |
|                               | Fluid line clamp(s) is closed             | Open clamp(s)                                                        |
|                               | Air bladder or internal air leak          | Service required                                                     |
| Air Bladder not deflate       | Control switch is not "OFF"               | Flip control switch to "DEPRESSURIZE"                                |
|                               | Rapid deflation Switch is "OFF"           | Flip rapid deflation switch to "ON"                                  |
|                               | Switch damaged                            | Service required                                                     |
|                               | Internal problem                          | Service required                                                     |
| Leakage of solution           | Bag not spiked securely                   | Secure spike in fluid bag                                            |
| Unable to adjust the pressure | Incorrect pressure at the air source      | Adjust the pressure of air source to 60 psi                          |
| Unable to reduce pressure     | Damaged regulation knob or pressure gauge | Service required                                                     |
|                               | Unit is not "OFF"                         | Turn control switch to "DEPRESSURIZE"                                |

### ◆ ENVIRONMENTAL CONDITIONS FOR STORAGE:

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

### ◆ EXPLANATION OF SYMBOLS:

|                                                                                    |                                      |                                                                                     |                                       |                                                                                     |                                                                                   |
|------------------------------------------------------------------------------------|--------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|
|  | Not made with natural rubber latex   |  | Manufacturer                          |  | Consult instructions for use (or consult electronic instructions for use).        |
|  | Caution                              |  | Serial Number                         |  | Federal(USA) law restricts this device to sale by or on the order of a physician. |
|  | Temperature limit                    |  | Medical device                        |  | Authorized representative in the European Community / European Union              |
|  | CE Marking according to MDR 2017/745 |  | UKCA Marking according to UK MDR 2002 |                                                                                     |                                                                                   |

 **LAGIS<sup>®</sup> ENTERPRISE CO., LTD.**

No. 29, Gong 1st Rd., Dajia Dist., Taichung City 437, Taiwan  
 TEL: 886-4-26820767 FAX: 886-4-26820766  
<http://www.lagis.com.tw>

 **Obelis s. a.**

Boulevard Général Wahis 53 1030 Brussels, BELGIUM  
 TEL: +(32)2. 732.59.54 FAX: +(32)2.732.60.03  
 E-Mail: mail@obelis.net

**Emergo Consulting (UK) Limited**  
 c/o Cr360 - UL International Compass House, Vision  
 Park Histon Cambridge CB24 9BZ, United Kingdom

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