www.lagis.com.tw/eifu_en.php

STERILE EO

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

◆ DEVICE DESCRIPTION :

Lagis Electrosurgical Suction Irrigation System comes in a pistol grip handle with two types of tubing (SAP-H, SAP-HT) and three types of shafts (SIP-340-L, SIP-340-N, SIP-340-S). The handles and shafts are completely interchangeable, which offers a variety of configurations to meet the surgeon's needs. This device has been sterilized by ethylene oxide prior to use and for single use only.


Lagis Electrosurgical Suction Irrigation System is compatible with commonly used suction and irrigation equipment and monopolar electrosurgical unit (ESU) that comply with IEC 60601-1 standards.

PRODUCT MODEL

Model Numbers	Description	Wire Length(cm)
SAP-H	Disposable hand control pistol grip handle with 11 cm tubing	290
SAP-HT	Disposable hand control pistol grip handle with 280 cm tubing for dual solution bags	290
SIP-340-L	Disposable L-hook electrode, 5 mm shaft, 340 mm length	-
SIP-340-N	Disposable needle electrode, 5 mm shaft, 340 mm length	-
SIP-340-S	Disposable spatula electrode, 5 mm shaft, 340 mm length	-

ILLUSTRATION AND NOMENCLATURE :

- 1.Dissection / coagulation buttons
- 2.Suction / Irrigation Buttons
- 3.Button pin
- 4.Power lead
- 5.Handle
- 6.Separation Button
- 7.Suction / Irrigation tubing connectors
- 8.Suction tubing connector
- 9.Tubing spikes
- 10.Connect part
- 11.Rotation knob
- 12.Insulating tube
- 13.Electrode

◆ INDICATIONS/INTENDED PURPOSE :



1. Lagis Electrosurgical Suction Irrigation System has application in minimally invasive procedures to facilitate tissue dissection, coagulation, irrigation, and fluid evacuation through all appropriately sized trocar sleeves.
2. The expected clinical benefits: reduce the surgery risks, shorten the surgery time, minimize the wound, and reduce the bleeding situation.
3. Patient target group: Lagis Electrosurgical Suction Irrigation System is intended for adults.

◆ CONTRAINDICATIONS :

1. This device is not intended for contraceptive coagulation on fallopian tissue. It may be used to achieve hemostasis following transection of the fallopian tube.
2. This device is not intended for use when minimally invasive techniques are contraindicated.

◆ INSTRUCTIONS FOR USE :

1. Insert the connect part to the handle until hearing the "click" sound.
2. Retract the button pin located at the right side of the handle to unlock the buttons.
3. Connect the suction tubing connector to a fluid collecting canister.
4. **SAP-H:** Connect the irrigation tubing connector to an irrigation tubing male connector.
SAP-HT: Insert the spikes into saline solution bags.
5. Connect the power lead to the ESU. Please make sure the ESU is compatible with the power lead prior to use.
6. When removing or inserting the device from the trocar sleeve, ensure that the insulating tube is pushed completely forward to cover the electrode tip of the device. This will prevent the device or trocar sleeve from any damage.
7. Introduce the shaft through the appropriate sized trocar sleeve to perform surgical operations:
 - Electrode tip exposure**-By sliding the rotation knob forward or backward to control the exposure of the electrode tip.
 - Electrode tip rotation**-Rotate the rotation knob to turn the electrode tip.
 - Suction/Irrigation control**-Press the suction button on the right to apply suction. Press the irrigation button on the left to irrigate. Generally, the appropriate pressure for irrigation is approximately 150-300 mmHg.

- Dissection/Coagulation**-Press the dissection button (yellow side) to dissect. Press the coagulation button (blue side) to coagulate.
- 8. To remove the rotation knob from the handle, press the separation button and pull the rotation knob away from the handle.
- 9. Once the surgery is completed, disconnect the device from the saline solution bag and the fluid collecting canister, and disconnect the power lead from the ESU.
- 10.Dispose of the used device properly and safely by healthcare professionals in accordance with the national waste management regulation.

◆ WARNINGS AND PRECAUTIONS :

1. Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
2. This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contraindications, and complications prior to performing the procedures.
3. This device is packaged and sterilized for single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocess, or resterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
4. Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
5. The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
6. In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in the instructions for use.
7. Regarding operating conditions:
 - Electrosurgery should not be performed in the presence of flammable solutions or gasses, or in high-oxygen concentration environments.
 - Studies have shown that smoke generated during electrosurgery may be harmful to surgical personnel. These studies recommend the use of a surgical mask and adequate surgical smoke evacuator or other means to prevent the smoke from inhaled by the surgical personnel.
8. Regarding equipment compatibility:
 - The shafts (SIP-340-L, SIP-340-N, SIP-340-S) are only compatible with the handles (SAP-H, SAP-HT), and vice versa. Do not attempt to connect these shafts and handles to any other electrosurgical connectors or instruments, which may cause serious injury to patients or users.
 - If trocars or accessories from different manufacturers are used together with Lagis Electrosurgical Suction Irrigation System, verify compatibility prior to performing the procedure and ensure that electrical insulation or grounding system is not compromised.
9. Regarding use of the equipment:
 - Do not use power generators that may deliver more than 2600 Volts peak to the Lagis Electrosurgical Suction Irrigation System. The electrosurgical generator power level should be set as low as is necessary to achieve the desired electrosurgical effect. Excessive power levels may result in instruments malfunction, insulation failure, and injury or thermal damage to patients or users.
 - Always use the lowest possible output setting necessary to achieve the desired electrosurgical effect. In order to lessen the possibility of unintended burn injury, use the active electrode only for the minimum time necessary. Higher settings may result in spurious insulation failure alarms and/or insulation breakdown.
 - If failing to achieve the desired electrosurgical effect with normal power level settings:
 - Check the grounding pad ("return electrode") is properly affixed to the correct position.
 - Check all electrosurgical cables are firmly connected to the ESU.
 - Make sure the shaft is fully inserted into the handle. Replace with another shaft if necessary.
 - Make sure the patient grounding pad ("return electrode") is affixed to the patient according to proper grounding techniques.
 - Users should regularly inspect the device for any damage, in particular, the electrosurgical cables.
 - In order to reduce the risks of unintended burn hazards, activate electrode tip only when it is necessary.
 - If the device is not able to function normally, do not attempt to repair the device in any way and replace it with a new one.
 - The device MUST NOT be used under the following conditions:
 - If the electrosurgical cable is loose or wrongly connected, or the metal of the cable connected to the handle is visibly exposed. It may result in patient or user injury.
 - If the connection of the electrosurgical cable and the handle is poor. It may result in patient or user injury.
 - If there are scratches or surface defects on the device, or the insulation of the shaft shows cracks. It may result in parts loosen and probably dropping into the field of operation, or even patient or user injury.

- Do not attempt to repair the Lagis Electrosurgical Suction Irrigation System. The use of a defective instrument or an instrument that becomes defective during surgery may result in loose parts that could drop into the field of operation.
- Improper assembly or operation of the device may result in prolong surgical time and tissue or organ injuries from mechanical or thermal damage to the patient or operator.

◆ **ENVIRONMENTAL CONDITIONS FOR STORAGE :**

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

◆ **EXPLANATION OF SYMBOLS :**

	Temperature limit		Caution		Authorized representative in the European Community
	Do not re-sterilize		Manufacturer		Not made with natural rubber latex
	Do not re-use		Consult instructions for use		Do not use if package is damaged and consult instructions for use
	Date of manufacture		Single sterile barrier system		Sterilized using ethylene oxide
	Medical device		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.		CE Marking according to MDR 2017/745

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**LAGIS[®] ENTERPRISE CO., LTD.**