
BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

◆ **DEVICE DESCRIPTION :**


The Multiple Instrument Access Port consists of three components: top housing, attachment ring, and LAGIS Wound Retractor. The top housing consists of combinations of different number of 5mm, 5-12mm, and 15mm ports. The stopcock valves on the top housing are compatible with luer lock to provide gas insufflation and desufflation. The top housing can be removed during the procedure, allowing for extraction of large specimens. The wound retractor is composed of two flexible rings interconnected by means of a sleeve. The top flexible ring allows connection to the attachment ring and top housing to provide circumferential wound retraction.

1. Top housing

2. Attachment ring

3. LAGIS Wound Retractor

◆ **PRODUCT MODEL:**

TRL-0110R, TRL-0130R, TRL-0210R, TRL-0211R, TRL-0220R, TRL-0310R, TRL-0400R, TRL80-0040R, TRL80-0130R, TRL80-0220R

◆ **INDICATIONS/INTENDED PURPOSE :**

- The Multiple Instrument Access Port is designed for use during laparoscopic surgical procedures for maintaining pneumoperitoneum while allowing for insertion of multiple surgical instruments through a single incision into the body cavity.
- Expected clinical benefits: reduce postoperative pain, better cosmetic results, a faster recovery and low morbidity by using single port laparoscopic surgery.
- Patient target group: the Multiple Instrument Access Port is intended for adults who need laparoscopic procedures.

◆ **CONTRAINDICATIONS :**

The Multiple Instrument Access Port is not intended for use when minimally invasive surgery techniques are contraindicated.


◆ **INSTRUCTIONS FOR USE :**

Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised.

◆ **INSTALLATION PREPARATION**

Make an incision between 2-6 cm to facilitate insertion of different size of the wound retractor through abdominal wall and peritoneum. Standard Hasson technique is recommended.


◆ **INSTALLATION AND ASSEMBLY :**


1 Prepare the surgical site according to standard procedure, making sure the skin is clean and dry.


2 Place the template over the incision site, and mark an appropriate length incision line using a sterile marker.


3 Make an incision along the marked incision line.

5 Gently grasp the retaining "green ring" of the LAGIS Wound Retractor at the 10 o'clock and 2 o'clock position and roll inward. Repeat until ultimate retraction is obtained.


7 Attach the top housing onto the attachment ring.


9 To insufflate, attach a gas line to the stopcock on the top housing and ensure the stopcock valve is open.


4 Insert the "blue ring" of the LAGIS Wound Retractor into the operative site.

6 Attach the attachment ring to the retaining green ring to make the wound retractor circular and completely self-retaining.


8 Engage the lever lock to firmly connect the top housing and the Attachment ring.


10 Start the procedure.

◆ **DISASSEMBLY :**

1 Disassemble the lever lock.

2 Disassemble the top housing from the attachment ring.

3 Disassemble the attachment ring from retaining LAGIS Wound Retractor.

4 Retrieve the LAGIS Wound Retractor by simply removing the blue ring from the peritoneal cavity.

◆ **WARNINGS AND PRECAUTIONS :**

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- If endoscopic instruments or accessories from different manufacturers are used together in a procedure, verify compatibility prior to performing the procedure.
- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or resterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Do not use the retractor where the incision length is greater than 6 cm as loss of pneumoperitoneum may occur.
- Covering the retractor and incision with a surgical drape during removal with umbilical tape or suture will minimize dispersion of bodily fluids outside of the surgical space.
- If loss of pneumoperitoneum occurs: 1) verify condition of retractor still provides adequate wound protection, and 2) place damp gauze under the retractor ring, between the lower ring and the fascia, to stop the air flow.

- The flexible length retractor functions as a wound protector. Integrity of the flexible length retractor should be verified prior to specimen removal or tissue contact. If the integrity of the flexible length retractor is suspected or cannot be verified, wound protection should not be assumed, and specimens should be removed using a retrieval pouch.
- Use caution when introducing or removing instruments or prosthetic mesh through the retractor in order to prevent inadvertent damage. Special care should be used when inserting sharp or angled edged endoscopic instruments to prevent tearing.
- Do not pierce the flexible length retractor or the flexible rings with a needle or install a suture into or through the device to facilitate removal of the flexible length retractor. Piercing the sleeve or ring of the flexible length retractor will result in damage and potential loss of pneumoperitoneum during the procedure. Separation of the flexible ring from the retractor sleeve may also result from this action.
- Please remove the Multiple Instrument Access Port top housing while introducing the gauze. Do not introduce the gauze into the cavity from the 5 mm or 12 mm ports to prevent damage of valves of the ports.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in this instructions for use.
- Improper assembly or operation of the device may result in prolong surgical time, patient allergy, and tissue or organ injuries from mechanical damage to the patient or operator.

◆ **ENVIRONMENTAL CONDITIONS FOR STORAGE :**

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

◆ **EXPLANATION OF SYMBOLS :**

	Temperature limit		Caution		Authorized representative in the European Community/European Union
	Do not re-sterilize		Manufacturer		Not made with natural rubber latex
	Do not re-use		Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Date of manufacture		Single sterile barrier system		Sterilized using irradiation
	Medical device		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.		CE Marking according to MDR 2017/745

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