
BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

◆ **DEVICE DESCRIPTION :**

The Transanal Access Port consists of top housing, attachment ring, access channel, and introducer. The top housing consists of combinations of different number of 5mm, 5-12mm, and 15mm ports. The stopcock valves on the top housing are compatible with luer lock to provide gas insufflation and desufflation. The top housing can be removed during the procedure, allowing for the extraction of large specimens. The access channel allows connection to the attachment ring and top housing to provide circumferential anal retraction. The introducer may be inserted to aid in placing the access channel into position.

For the device CRS-A01M without top housing, it is designed for use with the cap for wound retractor and trocars to maintain pneumoperitoneum.

1. Top housing

2. Attachment ring

3. Access channel

4. Introducer

TRANSANAL ACCESS PORT

Models	Top Housing	Attachment Ring	Access Channel	Introducer
CRS-60-0220M	O	O	O	O
CRS-A01M	X	O	O	O

◆ **INDICATIONS/INTENDED PURPOSE :**

- The Transanal Access Port is designed for use during transanal natural orifice transluminal endoscopic surgery for maintaining pneumoperitoneum while allowing for insertion of multiple surgical instruments through the anus into the body cavity.
- Expected clinical benefits: less postoperative pain, complications, wound infections, and hernia formation; improved cosmesis; lower short-term cost; and optimal midterm and long-term surgical and oncologic outcomes.
- Patient target group: the Transanal Access Port is intended for adults who need transanal minimally invasive surgery procedures.


◆ **CONTRAINDICATIONS :**

The Transanal Access Port is not intended for use when minimally invasive surgery techniques are contraindicated.


◆ **INSTRUCTIONS FOR USE :**

Prior to use, inspect the contents of the packaging and the packaging itself to ensure the integrity is not compromised.


INSTALLATION AND ASSEMBLY


Apply sterile lubrication to the distal end of the access channel.


Manually compress the distal end of the access channel in a folded form.


Use the forceps to clip the folded end of the access channel and place it into anus.


Make sure the distal end of the access channel is securely seated behind levator sling.


Apply the lubrication to the introducer to aid in placing the access channel into position.


Remove the introducer. Hold the access channel and suture the suture ports to secure.


Apply lubrication to the indentation of attachment ring and the access channel, and connect them.


Attach the top housing to the attachment ring.

Make sure the 2 lever locks are locked.


Attach insufflation tubing to one of the two stopcocks. Make sure the other stopcock valve is closed.


Start insufflation and the Transanal Access Platform is ready for use.


DISASSEMBLY


Disassemble the lever lock.


Disassemble the top housing from the attachment ring.

Disassemble the attachment ring from the access channel.


Remove sutures.


Gently pull out the access channel until it is removed from the anus.


CAP FOR WOUND RETRACTOR SERIES REPEAT STEP 1 TO 7


Use a needle like conical tip to create the first insertion within the boundary line.


According to the procedure, surgeons can insert other trocars by the same steps within the boundary line.

Attach the trocar-inserted Cap on the white ring.


Attach insufflation tubing to one of the two stopcocks. Make sure the other stopcock valve is closed.


Start the procedure.

DISASSEMBLY


Disassemble the trocar-inserted Cap from the attachment ring.


DISASSEMBLY

Disassemble the attachment ring from the access channel.

Remove sutures.

Gently pull out the access channel until it is removed from the anus.

◆ **WARNINGS AND PRECAUTIONS :**

- Do not use if the device or the sterile packaging is open or damaged, as this may cause an increase in the incidence of wound infections.
- If endoscopic instruments or accessories from different manufacturers are used together in a procedure, verify compatibility prior to performing the procedure.
- This device is packaged and sterilized for single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocess, or resterilization may create contamination, infection, cross-infection, or even compromise the integrity of the device and lead to device failure which in turn may result in patient injury, illness or death.
- This device should be operated only by surgeons with adequate training and familiarity with surgical techniques. Consult medical literature for further information regarding techniques, hazards, contra-indications, and complications prior to performing the procedures.
- Do not pierce the access channel with a needle or install a suture into or through the device to facilitate removal of the access channel. Piercing the access channel will result in damage and potential loss of pneumoperitoneum during the procedure.
- Please remove the top housing while introducing the gauze. Do not introduce the gauze into the cavity from the 5 mm or 12 mm ports to prevent damage of valves of the ports.
- Return the device and packaging to the distributor or implement safe disposal handling by the healthcare professionals if the sterile packaging is damaged or unintentionally opened before use.
- The used device which comes into contact with bodily fluids should be considered as biomedical waste and safely disposed of by healthcare professionals in compliance with the national waste management regulation to prevent patients, users, and other persons from biological contamination or infection.
- In the event of malfunction of the device, changes in performance that they may affect safety, or even any serious incident, please report to the manufacturer, the competent authority of the Member State and the authorized representative immediately according to the contact information provided on the label or in this instructions for use.
- Improper assembly or operation of the device may result in prolong surgical time, patient allergy, and tissue or organ injuries from mechanical damage to the patient or operator.

◆ **ENVIRONMENTAL CONDITIONS FOR STORAGE :**

Appropriate storage environment is a clean and dry area away from sunlight with a temperature range of 13~30°C (55.4~86°F).

◆ **EXPLANATION OF SYMBOLS :**

	Temperature limit		Caution		Authorized representative in the European Community/European Union
	Do not resterilize		Manufacturer		Not made with natural rubber latex
	Do not re-use		Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Date of manufacture		Single sterile barrier system		Sterilized using irradiation
	Medical device		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.		CE Marking according to MDR 2017/745

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LAGIS® ENTERPRISE CO., LTD.