

LAGIS Path™

Multiple Instrument Access Port

TRL series

For Effective Specimen Retrieval, Extracorporeal Anastomosis and Single Incision Surgeries

• Multi-Instrument Entry Channels

Designed to accommodate various diameters of laparoscopic instruments. Accommodates incisions from 1.5 cm to 3.5 cm and varying abdominal wall thicknesses.

• Detachable Buckle

Allows for easy and efficient specimen removal.

• Operating Platform

The flexible, transparent design enables precise instrument manipulation and clear observation while maintaining stable pneumoperitoneum.

• Integrated Port System Without Trocar Dependence

The integrated port eliminates the need for conventional trocar cannulas, and enabling a more streamlined and unobstructed surgical environment.



Code	Number of 5mm Ports	Number of 5~12mm Ports	Incision Length
TRL-0110R	1	1	2-4 cm
TRL-0210R	2	1	2-4 cm
TRL-0220R	2	2	2-4 cm
TRL-0310R	3	1	2-4 cm
TRL-0400R	4	0	2-4 cm

QTY/box : 1 pcs

Code	Number of 5mm Ports	Number of 5~12mm Ports	Incision Length
TRL80-0220R	2	2	3-6 cm
TRL80-0220R-N	2	2	3-6 cm

* TRL80-0220R-N: without attachment ring

QTY/box : 1 pcs



TRL80-0040R



TRL80-0040R-N

LAGIS LapShield™ Wound Retractor

360° Atraumatic Retraction & Wound Protection

- Distributes force evenly for optimal field of view
- Allows for maximum exposure with a minimum incision size
- Allows visualization of wound margins
- Soft materials and evenly dispersed tension promote patient post-operative comfort

Greater Exposure

- 360° of elastic tension maximizes working area providing significantly greater exposure

Easily Placed and Removed

- Inserted in seconds. Flexible ring slips gently into incision

Decreases Wound Infection Rate

- Helps protect wound from contamination and keeps wound edges moist



Code	Number of 5mm Ports	Number of 5~12mm Ports	Incision Length
TRL80-0130R	1	3	3-6 cm
TRL80-0040R	0	4	3-6 cm
TRL80-0040R-N	0	4	3-6 cm

* TRL80-0040R-N: without attachment ring

QTY/box : 1 pcs